

SLA 3D Printable Gas Diffusion Microextraction

Piassalonga, Gabriel Baroffaldi
`gabriel.baroffaldi@unesp.br`

de Lima Gomes, Paulo C. F.
`paulo.clairmont@unesp.br`

February 2025

1 Considerations

This Gas Diffusion Micro-extraction (GDME) device was projected to be printed using the SLA 3D printer type. An FDM or SLA 3D printer can print the acceptor solution holder.

Use instructions:

- Place the sample inside the printed GDME;
- Add the hydrophobic membrane on top of the GDME device;
- Place the acceptor solution holder above the hydrophobic membrane;
- Add the acceptor solution.

2 Dimensions

Figure 1: Dimensions of the 3D printed GDME device.




Figure 2: Back dimensions of the 3D printed GDME device.




Figure 3: Inside view of the GDME device.





Figure 4: Dimensions of the acceptor solution holder.



3 License

SLA 3D Printable Gas Diffusion Microextraction © 2025 by Gabriel Barofaldi Piassalonga, Paulo Clairmont Feitosa de Lima Gomes, São Paulo State University is licensed under Creative Commons Attribution-ShareAlike 4.0 International. To view a copy of this license, visit <https://creativecommons.org/licenses/by-sa/4.0/>